


knotilus

Suture Knot-Tying Instrument



Jessica Xu, Michael Gee,
Kent Shibata, Desmond Kuan

16-879 Medical Robotics
December 15, 2022


knotilus Proof-of-Concept

3 Electronics & Control

Looping 2

4 Housing & Integration

1 Needle Drive



KEY RESULTS


- ✓ Novel looping mechanism
- ✓ Clinically relevant knot quality
- ✓ Close to realistic needle size

3 Electronics & Control


4 Housing & Integration

1 Needle Drive


Looping 2



Components of Needle Drive Linkage System




Close-up: Needle Driver and Hook





Torsional spring connecting driver and hook for snapping motion of hook



Close-up: Needle Driver Track





Dual track allows needle driver to be normal to the needle surface at all times



Close-up: Needle Sleeve and Suture Thread Gap

Gap between base and needle sleeve
allows suture thread to not get tangled




3 Electronics & Control



4 Housing & Integration

1 Needle Drive



Looping 2



Hooks catch thread to make loops



Hook covers with optimal compliance




Too soft




Too hard




Just right



Hook and lip passively tilt down to release loops



Controlled suture release inspired by sewing machines




3 Electronics & Control

4 Housing & Integration

1 Needle Drive

Looping 2



Controls: Actuation Logic Flow

Looper forms loops
Servo 1



Send looper down
Linear Actuator




Drive needle through loops
Servo 2




Send looper up
Linear Actuator

On/Off Switch



'Next Action' Button

LED lights up when action completed and will not proceed until white button pressed.



Controls: Actuation Logic Flow

Looper forms loops
Servo 1



Send looper down
Linear Actuator




Drive needle through loops
Servo 2




Send looper up
Linear Actuator

Emergency Stop Button

Halt and retract needle if anything goes wrong.



Hardware: Circuit Diagram



LED

Digital OUTPUT, 5V

330Ω Resistor

Momentary Push Buttons
Continue, Emergency Stop



'Main circuit' 5V DC

Ground: Analog INPUT

On/Off Switch

In series with 'Main circuit'
5V DC

Hardware: Circuit Diagram



L289N Motor Driver

‘Main circuit’ 5V DC



- Power motor driver

‘Additional circuit’ 9V DC

- Power linear actuator

N3, N4: Digital INPUT

Hardware: Circuit Diagram





Servo Motors

'Main circuit' 5V DC

Feedback: Digital INPUT

Driver: Digital OUTPUT

Hardware: Circuit Diagram



Power Consumption

5V, Max Draw: 0.032A


Ideal for miniaturization and untethered handheld devices.

3 Electronics & Control

4 Housing & Integration


1 Needle Drive

Looping 2




Vertical looper movement to enter/exit needle path

Laser Cut Mount Board





Lateral Alignment



Retracted:
Waiting,
Looping

Extended:
In needle
path

Linear actuator



Back

Side view alignment

Laser Cut Mount Board

